lightforce®

LIGHTFORCE® Fxi

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




8.2. Laserproduktmærkat

		
		
		
	LiteCure LLC 101 Lukens Drive, Suite A New Castle, DE 19720 USA Tel: 302.709.0408 Pat. www.litecure.com/patents	 EMERGO EUROPE Prinsessegracht 20 2514 AP The Hague The Netherlands Phone: +31.70.345.8570
		
	This laser product complies with FDA performance standards for laser products except for conformance with IEC 60825-1 Ed. 3 and IEC 60601-2-22 Ed 3.1 as described in Laser Notice No. 36, dated May 8, 2019.	

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
1. Kundekontonummer: _____
2. Model (REF): _____
3. Serienummer (SN): _____

INSTRUCTIONS FOR USE

Product Family: LightForce® Laser Systems

Models:

LTC-1510-L-B6

LightForce® FXi

Read the entire User Manual before attempting to operate the device. The most current electronic version of this document can be found at the following URL:

<https://www.lightforcemedical.com/therapy-laser-user-manuals/>

IFU ID: RSR-000886-000 (5)

Effective Date: 2021-Nov-30

©2021 by LiteCure, LLC. All rights reserved.

No part of this document may be copied or reproduced in any form, including an electronic retrieval system or be made available on the Internet, a public network, or otherwise, without the prior written permission of LiteCure, LLC.

Table of Contents

1.	Safety	3	4.4.1.	Pigment Considerations.....	11
1.1.	Terminology	3	4.4.2.	Dose	11
1.2.	Laser Safety	3	4.4.3.	Setting Treatment Expectations	11
1.3.	Laser Treatment Controlled Area (LTCA)	3	4.4.4.	Treatment Technique.....	12
1.4.	Fire and Explosion Hazards	3	5.	Maintenance	12
1.5.	Safety Features.....	3	5.1.	Cleaning	12
1.5.1.	Emergency Stop Button	3	5.1.1.	Cleaning the Console	12
1.5.2.	Password Entry	3	5.1.2.	Cleaning the Large Massage Ball.....	12
1.5.3.	Remote Interlock Connector.....	3	5.2.	Before You Call—Troubleshooting.....	13
1.5.4.	Laser Emission Indicators	3	5.3.	Calibration Check Procedure	13
1.5.5.	Unattended Protection 100 Seconds.....	4	6.	Quality and Regulatory	14
1.5.6.	Internal Laser Energy Monitor	4	6.1.	Quality Systems	14
1.6.	Laser Eye Protection.....	4	6.2.	Device Classification	14
2.	System Set-Up	4	6.3.	Declaration of Conformity (EU)	14
2.1.	Location	4	7.	Specifications and Environmental Conditions	14
2.2.	Unpacking and Initial Set-up	4	7.1.	System Specifications (LightForce FXi)	14
2.3.	Battery Use and Charging.....	5	7.2.	Spare Parts	15
2.4.	Replacing the Battery.....	5	7.3.	Electromagnetic Compatibility (EMC) Tables.....	15
2.5.	Connecting the Beam Delivery Assembly.....	5	7.3.1.	Emissions Test Summary.....	15
3.	Device Description	6	7.3.2.	Immunity Test Summary.....	16
3.1.	Indications for Use	6	7.3.3.	Recommended Separation Distances	16
3.2.	Contraindications	6	7.3.4.	Guidance and Manufacturer's Declaration – Immunity – For Non-Life-Supporting Equipment	16
3.3.	Maximum Permissible Exposure (MPE)	6	7.4.	FCC	17
3.4.	Laser Console Diagrams	6	7.5.	ISED	17
3.5.	Empower™ Delivery System	7	7.6.	Radio Equipment Directive (RED).....	17
3.5.1.	Quick Disconnect Treatment Heads.....	7	7.7.	Disposal.....	17
3.5.2.	Treatment Head Attachments	7	8.	Labeling	17
3.6.	Fiber Management.....	8	8.1.	Explanation of Symbols.....	17
3.7.	Laser Safety Eyewear.....	8	8.2.	Laser Product Label	18
3.8.	Medical Grade Power Cord.....	8	8.3.	Laser Warnings	19
4.	Operation.....	8	8.4.	Warranty Seal.....	19
4.1.	Safety Precautions	8	8.5.	Unique Device Identification (UDI) Label	19
4.2.	Patient Preparation	9	9.	Contact Information	20
4.3.	Device Operation	9			
4.3.1.	Power On	9			
4.3.2.	Main Menu.....	9			
4.3.3.	Setup	9			
4.3.3.1.	Standard Settings.....	10			
4.3.3.2.	Language Setting	10			
4.3.4.	Protocols	10			
4.3.5.	Operation.....	10			
4.3.6.	Laser Emission.....	11			
4.4.	Laser Treatment Procedure	11			

1. Safety

1.1. Terminology

Various Warnings, Cautions, Recommendations, and Notes are presented throughout this manual and will be presented in the following manner:

Warning: Specific or potential danger. If ignored or compromised, the situation could result in serious injury.

Caution: Possible problem with the device associated with its use or misuse. Problems include, but are not limited to, device malfunction, device failure, or damage to the device.

Recommendation: Offers guidance for the optimal application and usage of the device.

Note: Describes the conditions or exceptions that may apply to the subject matter presented.

1.2. Laser Safety

Each facility and user of the laser device should follow all local laser safety requirements.

Federal, state, and local regulations as well as national and international standards should be considered when developing laser safety programs.

Each person using the laser system must have completed basic laser safety orientation and proficiency training prior to operating the laser system.

Recommendation: Designate at least one person at each facility that uses this laser system as the Laser Safety Officer, responsible for providing training on all operating and safety procedures and for monitoring and overseeing control of all laser hazards.

1.3. Laser Treatment Controlled Area (LTCA)

Select a secure and well-ventilated location in which to install and operate the laser system. The location should be "light-tight" (not have windows or viewing ports) and should have a single lockable door access. The access door should be closed and locked from the inside throughout the entire treatment with the laser system.

If the door is equipped with an interlock switch, the door can be connected to the remote interlock on the laser system – See Section 1.5.3. When the treatment door is opened during a treatment, the laser will stop emission, produce an error message, and return to **Standby** mode.

If you choose to set up your laser system in a multi-use area, a blocking barrier, screen, or curtain, certified for blocking the laser beam, must be used to create a smaller area inside a larger room. The barrier should be set up to simulate a "light-tight" room with single entry access as described above.

Warning: Always post a laser safety sign on the outside of the treatment room door or single access point when the laser is in use.

1.4. Fire and Explosion Hazards

Warning: Do not operate this laser system in areas that have explosion hazards such as flammable materials, gases or substances. A fire or explosion could occur. The laser beam can ignite most nonmetallic materials. A fire extinguisher should be readily available.

1.5. Safety Features

The laser system has multiple safety features. All individuals who use this laser system should be trained and made aware of the following features:

1.5.1. Emergency Stop Button

The button (**Figure 1**) is located on the front of the laser system and is illuminated when the machine is on. Push the button to turn off all electrical power to the device's microprocessor and stop the laser-emitting components.

Figure 1

1.5.2. Password Entry

The laser system is equipped with password entry as the key control, which is required after powering on the device, to prevent unauthorized system access.

1.5.3. Remote Interlock Connector

The laser system is equipped with a remote interlock connector (**Figure 2**) located on the bottom of the device. The remote interlock connector stops laser emission when the plug is removed or when the terminals of the connector are open.


Figure 2

The laser system is shipped with a mating connector (as shown in **Figure 2**). The mating connector **MUST** be inserted (or remote interlock circuit with closed door) for the laser system to emit.

1.5.4. Laser Emission Indicators

The state of the laser system is shown on the touch screen and front bezel and indicated through audible alerts and the aiming beam.

When the laser system is emitting, an EMISSION! indicator (**Figure 3**) will appear on the touch screen and flash. See Section 4.3.6 for a full description of the transitions of this indicator based on the laser mode (i.e. Standby, Ready, or Emission).

Figure 3

The amber LED (**Figure 4**) next to the power button on the front bezel illuminates to indicate the state of the laser. When off, the laser is considered to be in a safe state - **Standby**. While transitioning from **Standby** to **Ready** mode, the LED blinks. When in **Ready** mode, it is constantly illuminated. While emitting, it blinks.

Figure 4

In addition, there is a red aiming beam that emits from the hand piece while in **Ready** and **Emission** modes that provides a visible indication of laser state and the treatment location. This aiming beam flashes during the transition from **Standby** to **Ready** mode.

Lastly, there is an audible indication during **Emission** that is a steady tone, or continuous beep, depending on the system settings. The audible alert sounds six times during the six second transition from **Standby** to **Ready** mode.

1.5.5. Unattended Protection 100 Seconds

If the laser system is left in **Ready** mode without enabling laser emission, the system will automatically return to **Standby** mode.

1.5.6. Internal Laser Energy Monitor

To prevent over-exposure hazards due to excessive laser power, the laser system continuously monitors the output current of the diodes to ensure they are within the factory calibration settings. If it is detected that the electrical current levels exceed the upper or lower limits that would result in power levels outside of the $\pm 20\%$ thresholds, the laser will stop emission and return to **Standby** mode. To clear the error message and reset the laser system, the operator must press the **EXIT** button on the touch screen.

1.6. Laser Eye Protection

ANYONE in the Laser Treatment Controlled Area (LTCA) during laser treatment **MUST** wear the laser safety eyewear that is provided with the Laser System. The eyewear has an optical density (OD) rating >5.0 for 808 nm **AND** 980 nm wavelengths, as identified on the lenses (**Figure 5**). Additional protective eyewear can be purchased from the company webstore.

2. System Set-Up

Caution: ALWAYS turn the laser system **OFF** and unplug the power cord from the wall outlet before lifting, moving, or relocating the laser system.

2.1. Location

Select a secure and well-ventilated location in which to install and operate the laser system (see Section 1.3 for description of Laser Treatment Controlled Area).

- The device should be within 1.8 meters (6 feet) of an available 100-240 V electrical outlet.
- Select a flat hard surface that can adequately support the laser system.
- Ensure adequate airflow around the laser system. The laser system is air-cooled and designed for use in a well-ventilated location. There must be a minimum 10 cm (4 in) clearance around the back of the device.
- An appropriate fire extinguisher should be readily available.

Warning: Do not operate the laser system in areas with explosion hazards such as flammable materials or gases.

2.2. Unpacking and Initial Set-up

- Step 1.** Carefully unpack the laser system and its accessories from the shipping carton.
- Step 2.** Inspect the laser system carefully for any damage, such as cracks, dents or bent parts.
- Step 3.** If items are missing or damaged, call Customer Care. Also, please notify the carrier if the damage appears to have occurred during shipping and retain packaging for evidence.

Caution: If the Warranty Seal Label (**Figure 6**) is not found on the bottom of the device or is broken, **DO NOT** operate the device. Call Customer Care.

Figure 6

- Step 4.** Fully charge the battery before using this laser system. When connecting the external power supply to the laser system, make sure that the cord is correctly aligned. With the flat side of the cord facing up (**Figure 7**), plug the external power supply into the DC input on the back of the laser system (**Figure 8**).

Figure 8

Figure 7

Caution: Do not force the external power supply cord into the DC input in the incorrect orientation (e.g., upside down). This will damage the pins and the unit will not charge properly.

- Step 5.** Plug the AC power cord into the power supply.
- Step 6.** Once the power connection is secure, plug the AC power cord into a 100-240 V grounded outlet. A green light on the power supply should illuminate (signaling active source power connection). If the laser system is off, the battery will charge provided the laser system is correctly plugged into the power supply and the wall outlet. The charging time is approximately 3 hours.

Recommendation: Charge the battery before use. The laser system ships with the battery only partially charged.

2.3. Battery Use and Charging

Note: When the external power supply is connected to the laser system, the battery icon will blink green to indicate it is charging the battery. To determine the battery charge level, the power supply must be unplugged from the rear of the device.

Warning:

- Do not modify the battery or attempt to insert foreign objects into it.
- Do not short circuit the battery, keep away from metallic items.
- Do not immerse the battery in water or get the battery wet.

Caution:

- Use of this battery in other devices could result in damage to the device or battery.
- The temperature range for charging this battery is 10 °C to 32 °C (50 °F to 90 °F). Do not charge the battery outside of this range.
- Do not store the battery at temperatures above 70 °C (158 °F), such as inside a car on a hot day or in direct sunlight; this may damage the battery.
- Avoid dropping the battery. If you suspect the battery or device is damaged, contact Customer Care.
- Do not disassemble, open, crush, bend, deform, puncture, shred or burn the battery.

When the laser system is turned on, there is a battery icon in the upper right corner of the display screen. The battery icon will indicate the amount of charge in the battery when the laser system is not connected to AC power (mains/wall outlet). The icon will also change color from green to amber to red as the charge is depleted. When the battery icon is red and blinking, the battery needs to be recharged.

2.4. Replacing the Battery

Caution: Replace the battery only with an authorized replacement battery. These batteries have been tested in accordance with applicable safety standards.

Warning: Discarded batteries may create a safety hazard. Before disposal, cover the battery contacts with electrical insulating tape.

If the battery is not charging as expected, contact Customer Care to go through troubleshooting to determine if a replacement battery is needed.

If the battery needs to be replaced, follow the steps below for removal and replacement of the battery.

- Step 1.** Turn the laser over and place it on a flat surface. Remove the battery cover by turning the finger locks. (**Figure 9**)
- Step 2.** To remove the old battery, grab the removal strap and pull straight up. (**Figure 10**)
- Step 3.** With the arrow on the battery pointing to the front of the unit (**Figure 11**), insert the new battery into the battery compartment.


Figure 11


Figure 10

Figure 9

- Step 4.** With even pressure, fully seat the battery into the battery compartment.
- Step 5.** Replace the battery cover and secure the finger locks. The battery cover will fit tightly to hold the battery in place.

2.5. Connecting the Beam Delivery Assembly

Caution:

- Do not remove the hand piece fiber cable from the emission port once it has been secured (unless the optical fiber cable needs to be replaced).
- Repeated insertion and removal of the hand piece fiber cable will increase the chance of emission port or fiber tip contamination. If the emission port or the fiber tip is contaminated, the optical fiber can be damaged during laser emission.
- If the fiber must be removed, use the dust cap to help prevent dust and debris from contaminating the emission port. Do not leave the port unprotected.

Note: The laser system will be delivered with the optical fiber already connected. If it needs to be replacement, then follow Steps 1-6 to disconnect the old and connect the new handpiece to the console.

- Step 1.** Unscrew metal fiber connector until loose and gently pull fiber away from unit **Figure 12**.
- Step 2.** After removing fiber, unplug finger switch cord by pulling plug away from unit.
- Step 3.** Remove the new assembly from the protective foam packaging.

Figure 12

- Step 4.** Secure the finger switch connector by carefully pushing it into the finger switch port (**Figure 13**).
- Step 5.** Remove the protective dust cap from the end of the new optical fiber cable.
- Step 6.** Open the hinged door covering the emission port and insert the fiber tip straight into the emission port (**Figure 14**). Twist the metal connector clockwise onto the threaded emission port until tight.

Figure 14


Figure 13

3. Device Description

The laser systems include a therapeutic laser console and optical delivery system. The optical delivery system consists of a flexible optical fiber connected to a hand piece with interchangeable treatment heads. A custom software application allows the user to select from built-in treatment protocols that adjust dosage and wavelength based on the patient's characteristics and condition. The device also provides the ability to set custom treatment parameters including output power, treatment time, and pulse rate. Once the user has selected the desired treatment parameters, activation of laser emission requires that the operator follow a sequence of predetermined steps.

3.1. Indications for Use

810 nm and 980 nm wavelength:

The laser system emits energy in the visible and near-infrared spectrum to provide topical heating for the purpose of elevating tissue temperature for a temporary relief of minor muscle and joint pain and stiffness, minor arthritis pain, or muscle spasm; the temporary increase in local blood circulation; the temporary relaxation of muscle.

The device is meant for professional use only by licensed practitioners.

Caution: Federal, national, or local laws may restrict this device to sale by or on the order of a licensed practitioner licensed by the applicable authorities to use and/or prescribe use.

3.2. Contraindications

- Do not apply infrared light to abdominal or lumbosacral points in pregnant females.
- Do not apply infrared light to the epiphyseal lines in children.
- Do not apply infrared light to the thoracic region or over the pacemaker itself in patients with pacemakers.
- Do not apply infrared light over the thyroid gland, ovaries or testicles.
- Do not apply infrared light to patients who are taking drugs that have heat or light sensitive contraindications, such as but not limited to certain types of steroids.

3.3. Maximum Permissible Exposure (MPE)

Maximum Permissible Exposure (MPE) is the highest power density (in W/m²) from a light source that does not require clinical evidence to be considered safe. For the laser radiation from the laser system, the MPE for ocular and skin exposure is determined by the wavelength, beam divergence, and exposure time. For times between 10 and 3000 seconds, the ocular MPE is 418.9 W/m² and the skin MPE is 3027.1 W/m².

3.4. Laser Console Diagrams

1. Power Button
2. Emission Indicator LED
3. Emergency Stop Button
4. Touchscreen Display
5. Laser Warning Labels
6. Fiber Optic Cable
7. Strain Relief
8. Hand Piece Handle
9. Hand Piece Clip
10. Speaker
11. Fiber Manager Spool
12. USB Port
13. Remote Interlock
14. Foot Switch Port
15. DC Input (Power Supply Port)
16. Warranty Seal
17. Finger Switch Marking
18. Applied Part Marking
19. Optical Applicator Marking
20. Battery Cover
21. Product Label

3.5. Empower™ Delivery System

Caution: Do not bend the optical fiber too sharply: the minimum permissible bending radius is 40 mm. Mishandling of the fiber may lead to damage to the fiber or delivery system and/or harm to the patient or laser operator.

Warning:

- If any part of the Empower Delivery System overheats or produces smoke, immediately power off the Laser System, discontinue operation, and contact Customer Care for assistance.
- The use of accessories and/or optical fiber cables not specified for use with this Laser System may result in unsafe emissions or damage to the Laser System.

The Empower Delivery System is composed of a double-sheathed optic fiber and a hand piece with interchangeable treatment heads. The laser emission is controlled with a finger switch on the hand piece. During Emission, the LED indicator will illuminate. (Figure 15). The functionality of the finger switch is influenced by the “Override” setting in Setup. When the operator is administering a treatment with Override **Off**, the finger switch must be held down continuously to enable emission. When Override is **On**, the finger switch acts as an ON/OFF switch allowing the operator to press once and release to enable continuous emission.

Figure 15

Warning: The **Override On** setting for the finger switch should only be used by a qualified laser technician.

To end emission, the operator will then press and release again, which will return the laser to **Ready** mode. The laser system will remain in **Ready** mode until the operator resumes the treatment or selects **Standby** on the touchscreen to return to **Standby** mode.

Warning:

- **DO NOT** press the finger switch without first verifying the safe orientation and proper positioning of the distal end of the optical fiber and hand piece. Ensure compliance with all safety precautions.
- **DO NOT** operate the laser system without a treatment head (Figure 16). This can cause serious injury.
- **DO NOT** operate the laser system with an improperly attached or loose-fitting treatment head (Figure 17). This can cause serious injury.

Figure 17

Figure 16

Recommendation: As the aiming beam passes down the same delivery system as the working beam, it provides a good means of checking the integrity of the delivery system. If the aiming beam is not present at the distal end of the delivery system, its intensity is reduced or it looks diffused, this is a possible indication of a damaged or malfunctioning delivery system.

3.5.1. Quick Disconnect Treatment Heads

To attach a treatment head, locate the metal insert on the inside of the opening. Line up the slot within the metal insert with the bearing on the mating end of the hand piece (Figure 18). The hand piece and treatment heads are keyed so they will only fit and attach in this orientation.

Once they are lined up, push the attachment onto the end of the hand piece so it is fully seated. Then while holding the hand piece, rotate the attachment clockwise until it **clicks**. This indicates the attachment is fully engaged, locked in, and ready for use. (Figure 19)

Figure 18

Figure 19

3.5.2. Treatment Head Attachments

Beam divergence and the Nominal Ocular Hazard Distance (NOHD) for each of the treatment heads are listed in the below table. NOHD is defined in ANSI Z136.3 as “the distance along the axis of the unobstructed beam from the laser to the human eye beyond which the irradiance or radiant exposure during normal operation is not expected to exceed the appropriate Maximum Permissible Exposure (MPE).” At distances greater than the NOHD, the beam intensity becomes lower than the MPE.

Nominal Ocular Hazard Distance (NOHD) Table				
	Name	Image	Divergence (radians)	NOHD (meters)
OFF-CONTACT	Flat Window		0.815	0.6
	Small Cone		0.815	0.4
	Large Cone		0.815	0.4
ON-CONTACT	Small Massage Ball		0.226	1.1
	Large Massage Ball		0.057	4.9

3.6. Fiber Management

The side of the FXi console has an easy-to-use Fiber Manager Spool and Hand Piece Clip that allows the user to store the beam delivery assembly while not in use. When therapy is discontinued or between treatment cycles, the user may rest the hand piece in the Hand Piece Clip. For long term storage and during transport, the user may also wrap the fiber around the Fiber Manager Spool to avoid damage to the fiber.


Figure 20

When wrapping the hand piece fiber around the fiber manager spool, the first fiber wrap should be fed into the Fiber Slot so the plastic Bifurcator, which couples the finger switch and fiber optic cables, snaps into the slot. Five more wraps (for a total of six) should be made tightly around the spool to provide enough length for the hand piece to feed into the Hand Piece Clip (Figure 20).

3.7. Laser Safety Eyewear

Warning:

- **ALWAYS** wear the protective eyewear supplied with this laser system when operating the device.
- All personnel who are in the laser treatment controlled area **MUST** wear laser safety eyewear.
- **DO NOT** remove protective eyewear until the operator returns the laser device to Standby mode.
- If you suspect injury, such as direct exposure to the eyes, seek medical attention immediately.

Additional protective eyewear can be purchased from company webstore.

3.8. Medical Grade Power Cord

Caution: Do not use any power cord other than the one provided with the system.

If a replacement is needed, additional power cords can be purchased from the company webstore.

4. Operation

4.1. Safety Precautions

Anyone in the laser treatment controlled area during laser operation (including the operator, all assistants, and the patient) **MUST** wear the protective eyewear supplied with this laser system.

Recommendation: **DO NOT** allow any nonessential personnel into the laser treatment controlled area during any laser procedure.

Recommendation: **NEVER** allow untrained personnel to operate this laser system unless directly supervised by a properly trained and experienced individual.

Warning:

- The protective laser eyewear has an optical density (OD) rating > 5.0 for 810 nm AND 980 nm laser emission (see specification sheet with eyewear). **ONLY** wear the protective eyewear supplied with this laser system.
- All personnel in the laser treatment controlled area **MUST** wear laser safety eyewear.
- Replace any protective laser eyewear with eyewear from the manufacturer if there is any damage.
- **NEVER** look directly into the distal end of the optical fiber connected to an active laser device **WITH** or **WITHOUT** the appropriate laser-emission protective eyewear.
- **NEVER** aim the laser light directly into the eyes.
- **NEVER** direct the laser beam at anything other than the area to be treated. Direct or indirect eye contact with the output beam can cause serious damage, irreparable corneal and/or retinal damage, and possible blindness to one or both eyes.
- If you suspect injury, such as direct exposure to the eyes, seek medical attention immediately.
- **FAILURE TO COMPLY** with all safety instructions and warnings may expose everyone in the controlled laser treatment area to harmful levels of laser radiation and/or dangerous levels of electrical current.
- **FAILURE TO COMPLY** with the application techniques listed in the manual may lead to exposure to harmful levels of laser radiation.
- **AVOID THE USE** of flammable solvents, anesthetics, and oxidizing gases such as nitrous oxide (N₂O) and oxygen. The high temperatures produced during normal use of the laser system may ignite some materials, such as cotton or wool.
- After cleaning any treatment head with alcohol solution, **DO NOT** use until the alcohol solution completely evaporates. Doing so may cause the laser to ignite alcohol solutions or vapors.

Take the following steps to secure the treatment room, or the controlled area:

- Alert patients and personnel of laser safety precautions before they enter the laser treatment controlled area. Post a laser safety sign on the outside of the treatment room door when the laser is in use.
- Close the treatment room door during operation of the laser. The remote interlock connector can be connected to the treatment room door through an interlock circuit so laser emission stops when the treatment room door is opened.

4.2. Patient Preparation

Warning:

- Use carefully. May cause serious burns. Do not use over sensitive skin areas or in the presence of poor circulation. The unattended use of the device by children or incapacitated persons may be dangerous.
- **DO NOT** allow any reflective objects to fall into or obstruct the path of the laser while it is emitting. Direct or indirect eye contact with scattered laser light from any reflective surface can cause serious, irreparable corneal and/or retinal damage, and possible blindness to one or both eyes.
- **DO NOT** pretreat tissue (e.g., with a heating pad or wet compress) before laser treatment. Tissue to be treated must be dry and at normal temperature before laser treatment begins.

Before a treatment, do the following preparation:

- The operator, all assistants, and the patient must remove all reflective objects (rings, metal watchbands, and jewelry) prior to treatment.
- The area to be treated must be exposed. **DO NOT** treat through clothing or bandages.
- The treatment area should be clean and free of surface dirt or oils.
- Isopropyl alcohol solution (70% alcohol) may be used to clean all instrument surfaces that contact the patient (see Section 5.1). Remember to allow the alcohol solution to completely evaporate before initiating the Laser System.

4.3. Device Operation

4.3.1. Power On

After pressing the power button on the front of the console, there will be an audible beep and after initializing, the software displays the password screen (Figure 21). Press 1234 and the enter button.



Figure 21

4.3.2. Main Menu

The Main Menu screen can be accessed from other screens by pressing the home icon in the upper left corner.

4.3.3. Setup

In **Setup**, the software allows the operator to change the Aiming Beam, Unit Volume, and Override settings. From this screen, the user can Set the system Language.


4.3.3.1. Standard Settings

Aiming Beam – small pilot beam that provides the operator with a visual indication (red light) of the location of treatment beam. The aiming beam can be set to **Steady** or **Pulsed**.

Unit Volume – sets the sound that will occur during laser emission. The unit volume can be set to low, medium, or high with the option of a steady **Tone** or a **Beep**.

Override – sets how the operator controls the ON/OFF of laser emission. The Override is meant for advanced users who do not want to hold down the finger switch throughout the entire treatment. If the finger switch **VERRIDE** setting is ON, pressing the finger switch once will turn the laser emission on and laser emission will continue until the finger switch is depressed a second time (or until the time remaining for laser emission counts down to zero).

4.3.3.2. Language Setting

Set Language – enables the operator to set a language for all user interface screens. Pressing Set Language will take operator to the Language screen. On the Language Screen, select and highlight the desired language. Press "Save" button. The system will bring up a confirmation message for the language selected and will auto restart if operator selects "Yes". If the operator selects "No", the system will return to the Language selection screen.

4.3.4. Protocols

After selecting **Protocols** through the **Main Menu**, the Select Patient Characteristics screen (**Figure 22**) prompts user to select specific patient characteristics that can affect the absorption or penetration of light to the target tissue. The operator then selects the part of the body and condition to treat (**Figure 23**). The Influence Technology software will set the power, treatment time, pulse rate and recommended treatment heads to deliver the treatment. When using Protocol mode, the user cannot change the set time of the treatment.


Figure 23


Figure 22

4.3.5. Operation

When entering the **Operation** screen through the **Main Menu**, the operator is able to adjust the treatment time, power, pulse rate, and skin color. The laser will be in **Standby** mode and the **Operation** screen (**Figure 24**) will prompt the user to input all the treatment parameters. The aiming beam is not active in **Standby** mode.



Figure 24

The following adjustments and options are available from the **Operation** screen:

- **Time:** Press the time "+" or "-" button to adjust the laser treatment time. To scroll to the desired time value, maintain pressure on the desired button. The selected laser treatment time will be shown in the time window in minutes and seconds (Min:Sec).

- **Power:** Press the “+” or “-” buttons or rotate the power dial to adjust the laser power. The selected laser power will be shown in the power window in Watts as well as on the power dial. To scroll to the desired power value, maintain pressure on the dial and rotate in a clockwise direction. The user can also press on the squares around the perimeter of the dial to jump to higher or lower powers.
- **Pulse Rate:** Press the pulse rate “+” or “-” button to adjust the laser pulse rate. The selected laser frequency will be shown in the pulse rate window. To scroll to the desired frequency value, maintain pressure on the desired button. The laser frequency can be set in CW (continuous wave, not pulsed) or pulsed mode with frequencies of 2, 10, 20, 100, 200, 500, 1000, 2500, 5000, or 10,000 Hz. The duty cycle in the pulsed mode is 50%.

4.3.6. Laser Emission

When the operator has set the treatment parameters and is ready to start the treatment, the operator should press the **Standby** button to proceed to **Ready** mode.

The software has a **6 second delay (with 6 audible beeps)** to warn the user of the transition from **Standby** to **Ready** mode. During the transition, the **Standby** button changes to **Ready** and flashes. After 6 beeps, the laser is in **Ready** mode.

The purpose of **Ready** mode is to wait for the user to press the finger switch to start laser emission. In **Ready** mode, the laser will emit when the finger switch is activated. The laser system will return to **Standby** mode from **Ready** mode if emission is not enabled within 100 seconds.

Note: Pressing the finger switch during the 6 second delay will result in an **Operation Error**. The user will need to depress the finger switch prior to pressing exit, at which time the system will return to Standby.

 **Warning:** NEVER leave this device unattended in **Ready** mode.

In **Ready** mode:

- The **Ready** button is active and the aiming beam emits from the hand piece.
- The selected time, power, and pulse settings are still adjustable; however, the skin color is not adjustable.
- To start laser emission, press the finger switch. To stop emission, remove finger from the finger switch. The laser emission can be enabled again by pressing the finger switch until the remaining time is over.

Note: If the unit has the finger switch override ON, each time the finger switch is pressed, it will act as an on/off switch.

During emission, the word “**EMISSION!**” will replace the **Ready** button on the screen and the software does not allow the user to change any settings. In addition, the laser system will produce beeps, or a steady audible tone according to the audio settings chosen in **Setup**, and the aiming beam will continue to emit. Furthermore, the LED emission indicator on the front bezel will illuminate for the duration of the treatment.

4.4. Laser Treatment Procedure

This section contains information about the protocols and settings on this device. The following treatment recommendations are provided for guidance only. The practitioner is the person who determines the settings and protocol to use when treating each individual patient.

4.4.1. Pigment Considerations

This laser system will deliver a dual wavelength of 810 nm and 980 nm light which has been found to provide therapeutic relief of pain. This laser system allows the operator to choose skin color type for the patient. Since darker pigments will absorb 810 nm wavelength more than lighter pigments, the software will adjust the output so that it contains 100% 980 nm wavelength when the dark skin setting (Fitzpatrick V or VI) is selected.

Other circumstances when it may be appropriate to use this setting include:

- Highly pigmented areas, such as areas with a dark sun tan, or areas with a tattoo
- Areas with dense hair growth, such as the nape of the neck

4.4.2. Dose

The laser energy dose delivered during a treatment is measured in joules. This number is displayed on the screen during a treatment and is calculated by multiplying the average power (in Watts) by treatment time (in seconds). Laser treatment should be applied to an area using a scanning application.

The software protocols in the laser system are designed to deliver a dose appropriate for the treatment area and conditions input by the operator. It is recommended that treatment information (pre-set protocol and input settings or power, energy delivered, treatment time, and treatment area) is recorded for each treatment.

- In general, doses are delivered to the area of pain, the surrounding tissues, and along the nerve pathway for the specific area experiencing the pain.
- Treat muscles connecting to a painful joint in addition to tendon attachments.

4.4.3. Setting Treatment Expectations

Patients and clinicians need to understand the possibilities and limitations of the application of deep tissue laser therapy. Chronic pain may require several treatments before a beneficial effect is achieved. Deeper tissues may require multiple treatments before significant benefits are derived.

A typical treatment regimen may consist of three treatments spaced one day apart during the first week, two treatments the following week, and one treatment in week three. If no benefit is derived after the 6 to 10 treatments, then the patient’s options need to be reviewed considering the condition and response. Occasionally, a patient can experience a tingling sensation or soreness due to increased blood flow. In these instances, it is recommended to wait until the soreness subsides before continuing with additional treatments.

4.4.4. Treatment Technique

Warning:

- **DO NOT** allow any reflective objects to fall into or obstruct the path of the laser energy produced by this device. The operator, all assistants, and the patient must remove all reflective objects (such as rings, metal watchbands, and jewelry) prior to treatment with this device. Direct or indirect eye contact with scattered laser light from any reflective surface can cause serious, irreparable corneal and/or retinal damage, and possible blindness to one or both eyes.
- **DO NOT** pretreat tissue (e.g. with a heating pad or wet compress) before laser treatment. Tissue to be treated must be dry and at normal temperature before laser treatment begins.
- Use carefully. May cause serious burns. Do not use over sensitive skin area or in the presence of poor circulation.
- During treatment, it is important to get feedback from the patient concerning their comfort level. The laser provides a soothing warmth. Be cautious in cases where the patient has limited sensation or is not able to respond to increases in temperature and consider treating at a lower power setting. If the patient expresses any discomfort, the power output should be reduced.

- Treatment should be delivered to exposed skin. Do not treat through clothing or bandages. The treatment area should be clean and free of surface dirt or oils.
- Isopropyl alcohol solution (70% alcohol) may be used to clean all instrument surfaces that contact the patient (See Section 5.1 for cleaning instructions). Be sure to allow the solution to evaporate before initiating laser emission.
- For optimal delivery of the dose, keep the treatment head perpendicular to the skin surface. Continuously move the treatment head at a rate of approximately 3-10 cm/sec (1-3 in/sec).
- To ensure an even distribution of the dose to the treatment area, continuously move the treatment head over the area first moving continuously back-and-forth horizontally, and then moving back-and-forth vertically in a grid pattern (see **Figure 25**). Make sure to cover the entire area to be treated with this “painting motion”.
- Treatments can be applied in conjunction with stretching or gentle range of motion exercises.


Figure 25

5. Maintenance

Warning: **DO NOT** attempt to gain access to any internal component. Doing so may cause serious and/or irreversible injury. **THERE ARE NO USER-SERVICEABLE COMPONENTS** inside this laser device.

5.1. Cleaning

Warning:

- It is recommended to always turn off the laser system and unplug the power cord from the wall outlet before cleaning the system.
- Always use protective eyewear and gloves when cleaning and disinfecting any equipment.
- Clean the contact treatment heads (large massage ball and small massage ball) between treatments to prevent cross-contamination and infections.
- **DO NOT** use a treatment head immediately after cleaning with alcohol solution. Wait until the solution used in the cleaning procedure completely evaporates. Not allowing it to dry may cause the laser to ignite alcohol solutions or vapors.

Caution: **DO NOT** clean the lens inside the hand piece. Doing so may damage the lens during laser emission.

5.1.1. Cleaning the Console

Note: The laser system uses solid-state laser technology. It is important to keep the unit and accessories free from dust.

- Step 1.** Before cleaning the console, turn off the laser system and unplug the power cord from the wall outlet.
- Step 2.** The exterior surfaces on the unit and treatment heads can be cleaned by wiping with a soft, non-fibrous wipe (e.g., Kimwipes®) moistened with isopropyl alcohol solution.
 - Ensure the solution is only 70% alcohol. Solutions of more than 70% alcohol can cause product damage.
 - Avoid using soiled or gritty cloths.
- Step 3.** After cleaning, ensure that all cleansers have been removed and the parts are dry before use.

5.1.2. Cleaning the Large Massage Ball

- Step 1.** Remove the elastic band from the large massage ball treatment head.
- Step 2.** Unscrew the top portion of the large massage ball treatment head (**Figure 26**).
- Step 3.** The ball can now be removed for cleaning. Use a soft, non-fibrous wipe (e.g., Kimwipes®) moistened with 70% isopropyl alcohol solution. Do not wipe with a dry cloth, always dampen the cloth before wiping.

Figure 26

Recommendation: Be careful when unscrewing the top portion of the treatment head. It is recommended that this be done over a table to reduce the chance of dropping the ball directly on the floor.

- Step 4.** Ensure that all cleansers have evaporated and the parts are dry before re-assembly and use.
- Step 5.** Reassemble treatment head by replacing the massage ball and screwing the top portion securely on the treatment head.

Caution: Dropping the ball, may damage the ball and plastic housing. If significant damage has occurred, discontinue use. Ensure the ball is clean and free of debris that could cause damage over continued use.

If there are scratches on the massage ball:

There will be minor scratches on the massage ball resulting from normal use. These will not affect the efficacy of the laser treatment. **Figure 27** shows a massage ball treatment head with a ball that has some minor scratches.

Figure 28

Figure 27

Recommendation: If scratches in the large massage ball become significant (**Figure 28**), replace the massage ball treatment head.

5.2. Before You Call—Troubleshooting

If you are having a problem with your laser system, please check the list below of common conditions that can occur that you may be able to resolve without having to contact Customer Care.

The laser system will display error messages to caution the operator to fault conditions. If any of the following errors occur, the laser system will remain in **Standby** mode, or return to **Standby** mode, produce an audible warning and display the related message describing the fault condition. The following are examples of these fault conditions:

- Optical fiber is not inserted,
- Remote interlock is removed,
- Finger switch is depressed.

To proceed, correct the error and press the **exit** button.

The system continuously monitors the output current of the diodes to ensure they are within the factory calibration settings. If it is detected that the electrical current levels exceed the upper or lower limits that would result in power levels outside of the $\pm 20\%$ thresholds, the laser will stop emission and require system reset.

- When the laser pulse rate is CW (Continuous Wave), the laser output power is the same as the power setting on the display within operating tolerances.
- When the laser pulse rate is not CW (e.g. 10 Hz), the average output power is 50% of the power setting on the display, within operating tolerances.

The laser system will not enter or remain in **Emission** if the laser system is overheating. If temperatures exceed safety limits:

- There will be an audible warning.
- There will be an error message: Laser temperature is out of range.

Press the **exit** button and reset the system. Make sure that the ambient temperature of the laser system is within the operating temperature range of the device. If the problem still exists, contact Customer Care for further troubleshooting.

If the laser system will not turn on:

- Make sure that the power cord is correctly attached and plugged into a functioning power outlet.

5.3. Calibration Check Procedure

The instructions below provide steps for performing a calibration “check” that may be performed with certified, NIST traceable equipment. This procedure is recommended to be performed once a year but is not required unless there is evidence of significant damage or diminished output. If during the execution of this procedure, the stated requirements are not met, then factory maintenance and re-calibration may need to be performed on the device. Factory calibration must be performed by certified Service Personnel.

Caution: Calibration is a service procedure to be done only by certified Service Personnel. Adjustment to any internal components by anyone other than certified Service Personnel voids any existing manufacturer’s warranty on the Laser System.

Warning: Use of controls or adjustments or performance of procedures other than those specified herein may result in hazardous radiation exposure. Avoid eye or skin exposure to direct or scattered radiation when open.

Laser Power Calibration Check Instructions

Warning: ALWAYS wear laser safety goggles when performing the laser calibration check procedure and follow all laser safety guidelines.

Equipment Needed: Certified traceable power meter and detector (**Figure 29**) with appropriate wavelength and power measurement capabilities

- Step 1.** Turn off the laser.
- Step 2.** Inspect and attach optical fiber and power cord as directed. Make sure the optical fiber is clean and free of any dust, fluid or other contaminants.
- Step 3.** Turn on the laser system and enter Standby mode.
- Step 4.** Increase the power setting until it reaches the maximum wattage setting.
- Step 5.** Place the laser in Ready mode.
- Step 6.** Using the aiming beam, direct the distal end of the fiber into the active area of the power meter. Keep the fiber tip 2-3 cm away from the surface of the active area of the power detector.
- Step 7.** Activate the laser and record the value in watts from the power meter display. The laser power reading should be within $\pm 20\%$ of the power setting.
 - a. If the results are outside the 20% range, ensure that: all of the light from the fiber is entering the detector, the fiber is connected correctly, and the fiber is not damaged. Replace with a new fiber if necessary and repeat.
 - b. If the results are still outside the 20% range, discontinue this procedure and contact Customer Care for assistance.

Figure 29

6. Quality and Regulatory

6.1. Quality Systems

The FXi was designed, developed, and manufactured in accordance with ISO13485:2016 – Medical Devices – Quality management systems.

6.2. Device Classification

The FXi complies with 21 CFR Chapter 1, Sub chapter J, as administered by the Center for Devices and Radiological Health (CDRH) of the US Food and Drug Administration (FDA).

According to the applicable standards, the laser system is classified as follows:

- Class I Type B device per EN/IEC 60601-1
- Class 4 laser product according to IEC 60825-1

6.3. Declaration of Conformity (EU)

PRODUCT IDENTIFICATION	
Product name	Model/number
LightForce® Therapy Lasers	LTC-1510-L-B6 LightForce® FXi

MANUFACTURER		
Name of company	Address	Representative
LiteCure, LLC	101 Lukens Drive, Suite A New Castle, DE 19720 USA	Quality & Regulatory Director

AUTHORIZED REPRESENTATIVE		
Name of company	Address	Telephone/email
Emergo Europe	Prinsessegracht 20 2514 AP The Hague The Netherlands	+31.70.345.8570 - phone +31.70.346.7299 - fax europe@emergogroup.com

REGISTRATION INFORMATION		
Notified Body and ID #	Marking	CE certificate number
BSI Group, 2797		CE 542523

CONFORMITY ASSESSMENT		
Device classification	Route to compliance	Standards applied
Class IIa Rule 9	Annex II of MDD 93/42/EEC Council Directive	IEC 60601-1:2005+A1:2012 IEC 60601-1-2:2015 IEC 60601-1-6:2010+A1:2013 IEC 60601-2-22:2007+A1:2013 IEC 60825-1:2014

LiteCure, LLC declares that the above-mentioned product meets the fundamental requirements of the applied standards and the Medical Device Directive 93/42/EEC+2007/47/EEC as transposed in the national laws of the Member States.

7. Specifications and Environmental Conditions

7.1. System Specifications (LightForce FXi)

Device Family Name	LightForce Therapy Lasers - LTC	
Model Number	LTC-1510-L-B6	LightForce® FXi

Device Classifications	
USFDA Medical Device	Class II Medical Device per USFDA 21CFR 890.5500
Medical Electrical Equipment	Class I Type B device per EN/IEC 60601-1
Laser Product	Class 4 laser product per IEC 60825-1
User Interface	
Touch Screen Display	10" HD Display with Touch Screen
Emission Indicator	Audible Signal (50 to 75 dB), LED (front display)
Mechanical Specifications	
Dimensions	38 cm (L) x 28 cm (W) x 20 cm (H)
Weight	4 kg
Environmental Specifications	
Operation temperature	10°C to 35°C
Charging temperature	10°C to 32°C
Storage temperature	-20°C to 70°C
Humidity	≤ 80% RH Non-Condensing
Pressure	70-106 kPa
Cooling	Thermal Electrically Cooled with Forced Air

Optical Specifications		
Treatment Wavelengths	980 nm ± 20 nm & 810 nm ± 20 nm	
Skin Tone Adjustments	Light (I, II) or Medium (III, IV)	Dark (V, VI)
Maximum Output Power	15W	15W
Wavelength Ratio	80% of 980 nm + 20% of 810 nm	100% 980 nm
Aiming Beam Output	<4.0 mW	
Aiming Beam Wavelength	650 nm ± 20 nm	
Operating Modes	Continuous Wave (CW) or Pulsed	
Pulse Frequency	2, 10, 20, 100, 200, 500, 1000, 2500, 5000, 10000 Hz	
Pulse Duration	50% duty cycle at any given pulse frequency setting	
Laser Safety Goggles	OD5+ 808nm & 980nm	
Electrical Specifications		
Input Voltage	100-240 VAC; 50/60 Hz	
Input Current	≤ 400 VA	
Battery Pack	Li-Ion 9.9V, 5.0Ah	
Beam Delivery Assembly (BDA) Specifications		
BDA Model	New Empower Delivery System	
Hand piece Assembly (-B6)	Basic hand piece (BHP) with 600µm fiber core	
Empower Attachments	Flat Window-IQ, Small Cone-IQ, Small Massage Ball-IQ, Large Massage Ball-IQ, Large Cone-IQ (NOTE: "-IQ" indicates the attachment is the quick disconnect version with RFID tag)	

7.2. Spare Parts

Warning: Do not use any accessories and/or cables not specified or sold as replacement parts. Use of unauthorized accessories and/or optical fibers not specified or sold for use with this system may result in unsafe treatment emissions or damage to the laser system.

Item	Description	Image
Beam Delivery Assembly	Opaque protective double sheathed fiber with attached Empower™ hand piece and attachments	
Safety Goggles	808 nm and 980 nm OD 5+ Over glasses	
Safety Goggles	808 nm and 980 nm OD 5+ Regular	
Power Cord	Medical Grade Power Cord, 2.9 m	
External Power Supply	Medical Grade Power Supply	
Battery	49.5Whr, 9.9V, 5.0Ah (3S2P) Battery	

7.3. Electromagnetic Compatibility (EMC) Tables

This laser system has been tested to comply with the requirements of EN/IEC 60601-1-2 for electromagnetic compatibility (EMC). The following pages list the tests performed and corresponding test levels. The user, operator, installer, or assembler of this product is advised of the following:

- Medical electrical equipment needs special precautions regarding EMC and needs to be installed and put into service according to the information provided in this document.
- Portable and mobile RF communications equipment can affect medical electrical equipment.
- Use only the power cord provided with this product, or an approved alternate.

Warning: The use of accessories, transducers and cables other than those specified may result in increased emissions or decreased immunity of this product.

Warning: Equipment should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, this device should be observed to verify normal operation in the configuration in which it will be used.

All medical electronic devices must conform to the requirements of EN/IEC 60601-1-2. Precautions, adherences to the EMC guideline information provided, and verification of all medical devices in simultaneous operation are required to ensure the electromagnetic compatibility and coexistence of all other medical devices prior to a laser therapy treatment.

7.3.1. Emissions Test Summary

Specification	Frequency Range	Compliance Status
CISPR 11/EN 55011 Group 1, A Conducted Emissions	0.15 MHz – 30 MHz	PASS
CISPR 11/EN 55011 Group 1, A Radiated Emissions	30.0 MHz – 1,000 MHz	PASS

EN 61000-3-2 Power Line Harmonics	Up to the 40 th Harmonic	PASS
EN-61000-3-3 Power Line Flicker	Less than or equal to 4% Maximum Relative Voltage Change; Value of d(t) less than or equal to 3.3% for more than 500 ms	PASS

7.3.2. Immunity Test Summary

Specification	Minimum Test Level Required per EN 60601-1-2, 4 th Edition	Test Level Completed	Compliance Status
EN 61000-4-2 - Electrostatic Discharge Immunity	Air discharge up to ± 15 kV; Contact discharge up to ± 8 kV	Air discharge up to ± 15 kV; Contact discharge up to ± 8 kV	PASS
EN 61000-4-3 - RF Radiated Fields Immunity	Radiation field strength of 3V/m from 80 – 2700 MHz (80% AM, 1 kHz)	Radiation field strength of 3V/m from 80 – 2700 MHz (80% AM, 1 kHz)	PASS
EN 61000-4-4 -Electrical Fast Transient Immunity	Power line pulses of ± 2kV direct; I/O line pulses of ± 1kV	Power line pulses of ± 2kV direct; I/O line pulses of ± 1kV	PASS
EN 61000-4-5 -Lightning Surge Immunity	Power line surges of ± 2kV common, ± 1kV differential mode	Power line surges of ± 2kV common, ± 1kV differential mode	PASS
EN 61000-4-6 -RF Common Mode Immunity	150 kHz - 80 MHz at 3 Vrms 1 kHz 80% amplitude modulated	150 kHz - 80 MHz at 3 Vrms 1kHz 80% amplitude modulated	PASS
EN 61000-4-8 -Power Frequency Magnetic Field Immunity	Inductive loop at 50 Hz or 60Hz, to 30 amps (rms) per meter	Inductive loop at 50 Hz, to 30 amps (rms) per meter	PASS
EN 61000-4-11 - Voltage Dips and Short Interruptions	Voltage Dips of 100% for 0.5 cycles (0, 45, 90, 135, 180, 225, 270 and 315) 100% for 1 cycle 30% for 25 cycles 10% for 250 cycles (5s)	Voltage Dips of 100% for 0.5 cycles (0, 45, 90, 135, 180, 225, 270 and 315) 100% for 1 cycle 30% for 25 cycles 10% for 250 cycles (5s)	PASS

7.3.3. Recommended Separation Distances

The device is intended for use in the electromagnetic environment in which radiated disturbances are controlled. The customer or user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF Communications Equipment and the device as recommended below, according to the maximum output power of the communications equipment.

Below is the table of recommended separation distances between portable and mobile RF communications equipment and the device and systems that are not life-supporting.

Max Output Power (watts)	Separation (m) 150 kHz to 80 MHz $D=(3.5/V1)(\text{Sqrt } P)$	Separation (m) 80 to 800 MHz $D=(3.5/E1)(\text{Sqrt } P)$	Separation (m) 800 MHz to 2.5 GHz $D=(7/E1)(\text{Sqrt } P)$
0.01	0.1166	0.1166	0.2333
0.1	0.3689	0.3689	0.7378
1	1.1666	1.1666	2.3333
10	3.6893	3.6893	7.3786
100	11.666	11.666	23.333

7.3.4. Guidance and Manufacturer's Declaration – Immunity – For Non-Life-Supporting Equipment

The Laser System is intended for use in the electromagnetic environment specified below. The customer or user of the device should ensure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment – Guidance
Conducted RF IEC 61000- 4-6	3 Vrms 150 kHz to 80 MHz	(V1)Vrms	Portable and mobile communications equipment should be separated from the laser by no less than the distances calculated/listed below: $D=(3.5/V1)(\text{Sqrt } P)$ $D=(3.5/E1)(\text{Sqrt } P)$ 80 to 800 MHz $D=(7/E1)(\text{Sqrt } P)$ 800 MHz to 2.5 GHz where P is the max power in watts and D is the recommended separation distance in meters.
Radiated RF IEC 61000- 4-3	3 V/m 80 MHz to 2.5 GHz	(E1)V/m	Field strengths from fixed transmitters, as determined by an electromagnetic site survey, should be less than the compliance levels (V1 and E1). Interference may occur in the vicinity of equipment containing a transmitter.

7.4. FCC

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:

- (1) this device may not cause harmful interference, and
- (2) this device must accept any interference received, including interference that may cause undesired operation.

Caution: Any changes or modifications to this device not expressly approved by the manufacturer could void the user's authority to operate the equipment.

This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense.

This device complies with FCC RF radiation exposure limits set forth for an uncontrolled environment. This transmitter must not be co-located or operating in conjunction with any other antenna or transmitter. The equipment should be installed and operated with minimum distance of 20 cm between the radiator and your body.

7.5. ISED

This device contains license-exempt transmitter(s)/receiver(s) that comply with Innovation, Science and Economic Development Canada's license exempt RSS(s). Operation is subject to the following two conditions:

- (1) This device may not cause interference.
- (2) This device must accept any interference, including interference that may cause undesired operation of the device.

This Class B digital apparatus complies with Canadian ICES-003.

Radiation exposure statement:

This equipment complies with ISED radiation exposure limits set forth for an uncontrolled environment. This equipment should be installed and operated with greater than 20cm between the radiator & your body.

7.6. Radio Equipment Directive (RED)

The device is SAR exempt based on the operation of the device being under the SAR exemption requirements.

7.7. Disposal



If you plan to discontinue the use of this laser system, and intend to dispose of it or any of its parts or accessories, you must observe the applicable regional legal provisions for its disposal. You may also contact your local distributor, authorized service center, or Customer Care for the disposal options for this laser system.



8. Labeling

8.1. Explanation of Symbols


Symbol	Title	Description
	Manufacturer	Indicates the medical device manufacturer, as defined in EU Directives 90/386/EEC, 93/42/EEC and 98/79/EC.
	Date of manufacture	Indicates the date when the medical device was manufactured.
	Catalog number	Indicates the manufacturer's catalog, or model, number so that the medical device can be identified.
	Serial number	Indicates the manufacturer's serial number so that a specific medical device can be identified.
	Caution	Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.
	Type B Applied Part	Refers to the part of the medical device which comes into physical contact with the patient in order for the device to carry out its intended function.
	Do not dispose in unsorted municipal waste (WEEE)	Equipment must not be disposed of as unsorted municipal waste.
	Laser Warning	Warning label for class 2 and higher laser radiation
	Remote interlock connector	Identifies remote interlock connection port

	Foot Switch	Identifies connection port for foot switch
	Optical fiber applicator	Identifies connection port for hand piece fiber
	Emergency Laser Stop	Button used to terminate laser emission and shutdown the device in the event of an emergency.
	Fuse	To identify fuse boxes or their location. Note: Not user replaceable.
	Protective earth (ground)	To identify any terminal which is not intended for connection to an external conductor for protection against electrical shock in case of a fault, or the terminal of a protective earth (ground). Note: This is located inside the device.
	Follow Instructions	Indicates the need for the user to consult the instructions for use prior to operating the device.
	Power On/Off	Push/push power button
	USB Port	Connection for software updates and Instant Replay Backup / Restore
IOIOI	VGA Port	For Manufacturer Use ONLY
	Nemko-CCL Safety Mark with NRTL indicators	Indicates compliance with the Certification Body (Nemko-CCL) requirements regarding Electrical Safety (60601-1) in the US and Canada
	CE Marking of Conformity	Certification mark that indicates conformity with health, safety, and environmental protection standards for products sold within the European Economic Area.
	Prescription devices	CAUTION – Federal law restricts this device to sale by or on the order of a licensed practitioner licensed by the law of the State in which the practitioner practices to use or order the use of the device
	Authorized representative in the European Community	Indicates the Authorized representative (REP) in the European community (EC).

8.2. Laser Product Label

		
		
		
	LiteCure LLC 101 Lukens Drive, Suite A New Castle, DE 19720 USA Tel: 302.709.0408 Pat. www.litecure.com/patents	 EMERGO EUROPE Prinsesgracht 20 2514 AP The Hague The Netherlands Phone: +31.70.345.8570
		
	This laser product complies with FDA performance standards for laser products except for conformance with IEC 60825-1 Ed. 3 and IEC 60601-2-22 Ed 3.1 as described in Laser Notice No. 56, dated May 8, 2019.	

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8.4. Warranty Seal

Note: Warranty void if seal is broken.

8.5. Unique Device Identification (UDI) Label

Note: Example UDI shown for reference only. This format applies to medical, human use laser products only. For those models, the 2D, or QR, code represents the UDI text at the bottom of the label and each number provides information unique to the device as follows:

- (01) Indicates Device Identifier number, which is a Global Trade Identification Number (GTIN-14) that can be used to search for the device registration in the FDA GUDID Database.
- (11) Indicates Date of Manufacture in YYMMDD format.
- (21) Indicates Serial Number.

The linear barcodes are for manufacturer use only.


9. Contact Information

For immediate assistance, contact Customer Care directly.

If this laser system does not function as expected and/or if the laser system malfunctions, contact Customer Care using the contact information on the back cover.

Before you call, please have the following information so that Customer Care can provide you with the highest level of service:

1. Customer Account Number: _____
2. Model (REF) : _____
3. Serial Number (SN) : _____

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DJO Global
1a Guildford Business Park
Guildford, Surrey
GU2 8XG
UNITED KINGDOM

djoglobal.eu

© 2021 DJO · RSR-000886-000 · MKT00-11192-Rev A
Individual results may vary. Neither DJO Global, Inc. nor any of its subsidiaries dispense medical advice. The contents of this brochure do not constitute medical, legal, or any other type of professional advice. Rather, please consult your healthcare professional for information on the courses of treatment, if any, which may be appropriate for you. Deep tissue laser therapy™ is a medical procedure that is administered by a trained healthcare provider.

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U.S.A.